

**Item Name :** Water Bath

**Item Code :** 4263.2008.01

**Item Description :**

Double walled construction with inner chamber made of stainless and outer wall made of M.S. sheet Duly powder coated painted. The gap between outer & inner walls is filled with special grade glass wool insulation. The concentric rings can be removed to accommodate containers of various sizes. Temp. from ambient to 95°C ±2°C is to controlled by thermostat. Made out of heavy gauge S.S.304 with concentric rings. Works on 220/230 volts A.C. supply

**Salient Features:**

- Novel double walled structure effectively prevents vapor leakage.
- Intelligent temperature controller with PID program (Optional)
- Digital display and touch operating buttons (Optional)
- Function of timing and over temperature protection (Optional)
- Stainless Steel inner chamber



**Item Image :**

Size inside (Stainless Steel)	No. of Holes	Cap.(approx.)
300 x 250 x 100 mm	6 Holes of 75 mm dia.	8 ltrs.
355 x 405 x 100 mm	12 Holes of 75 mm dia.	15 ltrs.