

Cat No.

2067.2008.01



**GLACIER INDIA**  
**(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)**



**Item Name :** Litmus Paper (Red)

**Item Code :** 2067.2008.01

**Item Description :**

The red litmus will turn blue when it comes into contact with an alkaline solution above pH 8.3. Supplied in plastic snap cap container, about 100 paper strips per container. Litmus paper is NOT used to determine the pH level of a solution. If you need to find out the pH level, the correct test paper to use is the pH paper. • Pkt of 100 Strips

**Item Image :**



# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>