

**Item Name :** Slides Cabinets 1000 Slides

**Item Code :** 4259.2008.01

**Item Description :**

Wooden cabinet for keeping standard 75 x 25 mm glass microscope slides in vertical position. Each of the 10 drawers holds 100 slides in wooden groves, slides rest on a felt liner. Each drawer has a label holder and the outer door locks for security.

- Sturdy wooden construction • 1000 Slide Capacity • 10 drawers
- Stores slides vertically • Outer door locks



**Item Image :**

Cat. No.

4259.2008.01