

Item Name : Benchtop pH METER Digital

Item Code : 4250.2008.01

Item Description :

• Display 3½ LED • Range 0-14pH, 0±1999mV • Resulation 0.01pH, 1mV • Temperature Compensation Automatic • Power Requirement 220VA

Item Image :



Cat. No.

4250.2008.01

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>