

**Item Name :** Universal Hot Air Oven Memmert type

**Item Code :** 4245.2008.01

**Item Description :**

Standard double wall fabrication. Inner chamber made of richly anodized aluminum or highly polished stainless steel sheet. Exterior fabrication. The gap between the walls is filled with special grade glass-wool for heat insulation. Inner chamber is ribbed for placing shelves at convenient levels. Provided with two or three removable shelves. Special elegant door closing device has been provided. Nichrome wire heating elements - provided on three sides attain quick and uniform heating in range of 50°C to 250°C + 1°C controlled by capillary type thermostat. L-shaped thermometer is built-in type. Control panel is provided with selector switch of high or low rates of power, thermostat control knob and indicators for mains & thermostat. Supplied with cord and plug. The equipment is suitable to operate on 220 V AC 50 Hz single phase.

Chamber size in mms/inches) No. of shelves

300 x 300 x 300 (12" x 12" x 12" ) 2

Optional :- • Air Circulating fan.

- Timer 0-24 hours • Digital Temp. Indicator-Cum- Controller.
- Imported Thermostat "EGO"/ "JUMO" German with the accuracy of  $\pm 1^\circ\text{C}$ .
- Digital temp. Indicator.

**Item Image :**



Cat. No.

4245.2008.01

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>