

Cat No.

1167.2008.01



GLACIER INDIA
(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)



Item Name : Meker Burner

Item Code : 1167.2008.01

Item Description :

Meker burners are perfect for high-temperature work, and this one is capable of using natural gas. The nickel-plated brass tube features a non-removable air regulator that can be turned to adjust how much oxygen is feeding the flame, allowing you to adjust the height and intensity. The wire mesh across the mouth of the Meker burner allows for an evenly aerated flame of uniform intensity, due to the more intimate mix of air and gas.

Item Image :



44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>