

Item Name : Leaf Area Meter

Item Code : 4212.2008.01

Item Description :

Handheld leaf Area Meter is a portable instrument for the non-destructive measurement of leaf area. it can be adapted to have a conveyor for the continuous measurement of a large number of leaf samples. leaf Area Meter uses advanced laser technology to enable researchers a precise and convenient way to measure leaf area (or the areas of leaf-like objects).

Features

- Six measurement parameters: area, width, length, perimeter, shape factor, and aspect ratio
- Durable and lightweight • Improved consistency of measurements
- Improved handling of delicate leaves • SD Card for data storage
- Graphic display and rendering of leaf images for verification of a successful scan and auditing of data
- Built in GPS tagging to provide location data for each measurement
- Standard USB charging and data download • Ideal for field work and very easy to use
- Non-destructive measurements • No calibration required
- Root Length Measurement Attachment) is available
- Rechargeable battery
- Measures objects up to 150mm (254mm with conveyor) wide and 25mm thick of virtually unlimited length

Item Image :



Cat. No.

4212.2008.01

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>