

**Item Name :** Insect collecting net

**Item Code :** 4210.2008.01

**Item Description :**

- A) Insect Collecting Net With Aluminum pipe handle with fine rubber grip at end,length of handle 75 cm. 75 cm long net;fitted in ring through 5 cm wide cloth rim. Dia of bag at free end 25 cm. Ring dia 30 cm.
- B) Insect Collecting Net With Iron Handle : Like but with iron, chrome plated, pipe handle.
- C) Insect Collecting Net With Wooden Handle : Like but with wooden, spirit polish handle.
- D) Spare net bags for insect collecting nets.

**Item Image :**



Cat. No.

4210.2008.01