

**Item Name :** Incubator Stainless Steel

**Item Code :** 4209.2008.01

**Item Description :**

Standard double wall fabrication. Inner chamber made out of richly anodized aluminium or highly polished stainless steel sheet and outer insulation. Double wall door with double window for observation in the chamber. Temperature range 5°C above ambient to 80°C + 2°C control knob etc. The equipment is workable on 220 VAC 50 Hz single phase. Heating is achieved from 3 sides. Perforated removable shelves. Chamber size in mm. & Inches 300 x 300 x 300 (12" x 12" x 12")

**Item Image :**



Cat. No.

4209.2008.01