

Cat No.

2073.2008.01



**GLACIER INDIA**  
**(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)**



**Item Name :** Optical Spectrometer

**Item Code :** 2073.2008.01

**Item Description :**

Precision spectrometers have higher accuracy and trouble free workmanship, all mechanical and optical parts precisely manufactured and on table and movements is completely independent and has lock screw and slow movement knobs, so that movement can be read with the

- Collimator has rack & pinion focusing system and fitted with 7" focus Ach.
- Objective : Fine jaw adjustable slit
- Prism table rotatable to 360° and adjustable in height and have 3 leveling screws
- Spectrometer 6" scale with 1 minute verniers

Standard Accessories:

- Greeting holder • Prism holder • Reading magnifier
- Dustcover • Instruction manual • Plywood box

**Item Image :**



# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>