

Cat No.

2075.2008.01



GLACIER INDIA
(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)



Item Name : Hot Air Oven

Item Code : 2075.2008.01

Item Description :

(Three side heating elements) Standard double wall fabrication.

Inner chamber made of richly anodized aluminum or highly polished stainless steel sheet.

Exterior fabricated out of thick mild steel sheet duly finished in stove enamel off-white paint or paste mat finish colour combinations. The gap between the walls is filled with-special grade glass-wool for heat insulation. Inner chamber is ribbed for placing shelves at convenient levels. Provided with two or three removable shelves. Special elegant door closing device

has been provided. Nichrome wire heating elements -provided on three sides attain quick and

uniform heating in range of 50°C to 250°C + 1°C controlled by capillary type thermostat. Lshaped thermometer is built-in type. Control panel

low rates of power, thermostate control knob and indicators for mains & thermostat. Supplied

with cord and plug. The equipment is suitable to operate on 220 V AC 50 Hz single phase.

Chamber size in mms/inches) No. of shelves

350 x 350 x 350 (14" x 14" x 14") 2

Optional:- • Digital Temp. Indicator-Cum- Controller

• Air Circulating fan. • Timer 0-24 hours

• Digital Temp. Indicator-Cum- Controller • Digital temp. Indicator

• Imported Thermostat with the accuracy of $\pm 1^{\circ}\text{C}$.

Item Image :



A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>