

Cat No.

2077.2008.01



GLACIER INDIA
(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)



Item Name : Ph Meter Electrical Pen Type

Item Code : 2077.2008.01

Item Description :

A lightweight pocket sized pH tester providing fast and accurate reading on a large liquid crystal display.

- Range : 0 to 14pH • Resolution : 0.1pH
- Accuracy : +/- 0.1 pH • Weight : 85gm

Item Image :



44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>