

Cat No.

433.2008.01



**GLACIER INDIA**  
**(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)**



**Item Name :** Quickfit Apparatus

**Item Code :** 433.2008.01

**Item Description :**

This distillation apparatus combines a normally separate Liebig condenser, still head and receiver adapter into a one piece assembly.

- Comprises Liebig condenser, still head and receiver adapter in a single piece
- Angle and position pre-defined • Simple to set up with fewer clamps required
- Less grease contamination • No air leakage under vacuum
- No liquid hold-up at the joints

**Item Image :**



# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>