

Cat No.

434.2008.01



GLACIER INDIA
(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)



Item Name : Specific Gravity Bottle

Item Code : 434.2008.01

Item Description :

Specific Gravity Bottle with stopper of 50ml capacity, which is available in twopatterns:

- The ordinary capillary type specific gravity bottle with a neck of 6mm diameter.
- The wide mouthed capillary type specific gravity bottle with a neck of 25mm diameter.
- The ordinary specific gravity bottle is used for materials which remain absolutely fluid at 27°C while the wide mouth capillary type is used for materials which remain semisolid of high viscous at 27°C.

Item Image :



44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>