

Cat No.

361.2008.01



**GLACIER INDIA**  
**(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)**



**Item Name :** Thermometer

**Item Code :** 361.2008.01

**Item Description :**

These thermometers are ideal for use in laboratory, biology, chemistry, university, photographic, food/beverage, petrochemical, automotive, that has a dark scale contrast against a white background. Tech Specification

- Product Type: Glass Thermometer • Min temperature (0C): -20
- Max temperature (0C): 110 • Division (°C):1
- Length (mm):300 • Material: Glass
- Accuracy:  $\pm 1.5^\circ$ ,  $\pm 2^\circ$  above 100°C • Fill : Red Spirit

**Item Image :**



# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>