

Cat No.

2080.2008.01



GLACIER INDIA
(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)



Item Name : Spectrometer Kirchhoff-Bunsen Type

Item Code : 2080.2008.01

Item Description :

For qualitative observation and measurement of emission and absorption spectra. Table unit has fixed optical flint prism with removable cover. Includes swiveling observation telescope with movable ocular, scale tube with reference subdivision, slit tube with adjustable slit. Complete assembly supported on vertical pillar mounted on circular cast metal base.



Item Image :

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>