

Cat No.

2081.2008.01



GLACIER INDIA
(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)



Item Name : Standard Spectrometer

Item Code : 2081.2008.01

Item Description :

An economically priced instrument, which is nevertheless capable of much useful quantitative work. The main structural parts, including the telescope, collimator, and scale, are mounted on a common base with locking screw. A spring-loaded vernier scale attached to the telescope mount provides readings to 0.1° (6 minutes of arc). Collimator: 6mm long. Telescope: Mounted on moveable pillar with adjustment, locking screw and axis adjustment. Objective lens has 170mm.FL, 25mm diameter. Base: Mounted on a circular base with three leveling screws and has lined marked to assist placement of prism. Standard Accessories Supplied: 1 prism clamp for prisms up to 25mm diameter.
Specifications: • Scale: 170 mm diameter • Aperture: 25 x 25 mm



Item Image :

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>