

Cat No.

2083.2008.01



**GLACIER INDIA**  
**(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)**



**Item Name :** Sprayer

**Item Code :** 2083.2008.01

**Item Description :**

- Made from chemical resistant HDPE • Shape: cylindrical
- Lab bottles valve can be cleaned swiftly
- Used to dispense thin liquids
- Includes a plastic strainer, polypropylene trigger, intake and a valve



**Item Image :**

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>