

**Item Name :** Sockets Full Length Single

**Item Code :** 004.2008.01

**Item Description :**

**Item Image :**



Cat. No.	Socket Size	Approx. O.D. of tube mm.	Minimum Shank Length mm
004.2008.01	5/20	8	100
004.2008.02	7/25	10	100
004.2008.03	10/30	13	100
004.2008.04	12/32	16	100
004.2008.05	14/35	18	100
004.2008.06	19/38	22	100
004.2008.07	24/40	28	100
004.2008.08	29/42	32	100
004.2008.09	34/45	38	100
004.2008.10	45/50	50	100

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>