

**Item Name :** Expansion Adapters

**Item Code :** 014.2008.01

**Item Description :**

**Item Image :**



Cat. No.	Socket Size	Cone Size
014.2008.01	19/26	14/23
014.2008.02	24/29	14/23
014.2008.03	24/29	19/26
014.2008.03A	29/32	14/23
014.2008.04	29/32	19/26
014.2008.05	29/32	24/29
014.2008.06	34/35	19/26
014.2008.07	34/35	24/29
014.2008.08	34/35	29/32
014.2008.09	40/38	24/29
014.2008.10	45/40	29/32

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>