

**Item Name :** Multiple Adapters Two Necks Parallel

**Item Code :** 015.2008.01

**Item Description :**



**Item Image :**

| Cat. No.     | Socket Size | Cone Size |
|--------------|-------------|-----------|
| 015.2008.01  | 14/23       | 14/23     |
| 015.2008.01A | 14/23       | 19/26     |
| 015.2008.02  | 19/26       | 19/26     |
| 015.2008.03  | 19/26       | 24/29     |
| 015.2008.04  | 24/29       | 24/29     |
| 015.2008.05  | 19/26       | 34/35     |
| 015.2008.06  | 29/32       | 29/32     |