

Item Name : Recovery Bend Sloping

Item Code : 027.2008.01

Item Description :



Item Image :

| Cat. No. | Cone Size to fit Flask | Cone Size to fit Condenser |
|-------------|------------------------|----------------------------|
| 027.2008.01 | 14/23 | 14/23 |
| 027.2008.02 | 24/29 | 14/23 |
| 027.2008.03 | 19/26 | 19/26 |
| 027.2008.04 | 24/29 | 19/26 |
| 027.2008.05 | 29/32 | 19/26 |
| 027.2008.06 | 29/32 | 29/32 |
| 027.2008.07 | 24/29 | 24/29 |
| 027.2008.08 | 34/35 | 24/29 |
| 027.2008.09 | 29/32 | 24/29 |
| 027.2008.10 | 34/35 | 34/35 |

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>