

**Item Name :** Adapters Screwthread

**Item Code :** 040A.2008.01

**Item Description :**

Double-ended screwthread joints, useful for connecting to glass tubing etc. Supplied Complete with caps, rings and washers.

**Item Image :**



Cat. No.	Thread Size	Approx Length of limbs mm	Acceptable Dia. through screwcap mm
040A.2008.01	13	85	6.0-7.0
040A.2008.02	18	85	7.0-8.5
040A.2008.03	24	100	10.5-11.5
040A.2008.04	24	230	6.0-7.0
040A.2008.05	28	120	18.0-19.0