

Item Name : Flask Round Bottom Single Neck

Item Code : 058.2008.02

Item Description :

Clear glass, DIN 12348 Short

Note: If you need flasks in short & long height please mention separately

Item Image :



Cat. No.	Nominal capacity ml	Socket Size	Approx. overall height mm.	Dia mm
058.2008.02	5	14/23	65	27
058.2008.04	10	14/23	65	33
058.2008.05	25	14/23	85	41
058.2008.06	25	19/26	85	41
058.2008.07	25	24/29	85	41
058.2008.08	50	14/23	105	51
058.2008.09	50	19/26	105	51
058.2008.10	50	24/29	105	51
058.2008.10A	50	29/32	105	51
058.2008.11	100	14/23	115	64

058.2008.12	100	19/26	115	64
058.2008.13	100	24/29	115	64
058.2008.14	100	29/32	115	64
058.2008.16	150	19/26	120	74
058.2008.17	150	24/29	120	74
058.2008.18	150	29/32	120	74
058.2008.20	250	14/23	145	85
058.2008.21	250	19/26	145	85
058.2008.22	250	24/29	145	85
058.2008.23	250	29/32	145	85
058.2008.24	250	34/35	145	85
058.2008.25	500	19/26	175	105
058.2008.26	500	24/29	175	105
058.2008.27	500	29/32	175	105
058.2008.28	500	34/35	175	105
058.2008.29	1000	19/26	210	131
058.2008.30	1000	24/29	210	131
058.2008.31	1000	29/32	210	131
058.2008.32	1000	34/35	210	131
058.2008.33	2000	24/29	260	166
058.2008.34	2000	29/32	260	166
058.2008.35	2000	34/35	260	166
058.2008.36	5000	34/35	280	223
058.2008.37	5000	40/38	280	223
058.2008.38	10000	34/35	385	279
058.2008.39	10000	45/40	385	279

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>