

Item Name : Flasks Flat Bottom

Item Code : 059.2008.01

Item Description :

Single Neck DIN 12348 Medium

Item Image :



Cat. No.	Nominal capacity ml	Socket size	Approx. overall height mm.	Dia mm
059.2008.01	50	19/26	90	51
059.2008.02	50	24/29	90	51
059.2008.02A	50	29/3	90	51
059.2008.03	100	19/26	110	64
059.2008.04	100	24/29	110	50
059.2008.05	100	29/32	110	64
059.2008.07	150	19/26	115	74
059.2008.08	150	24/29	115	74
059.2008.09	150	29/32	115	74
059.2008.11	250	19/26	140	85

059.2008.12	250	24/29	140	85
059.2008.13	250	29/32	140	85
059.2008.14	250	34/35	140	85
059.2008.15	500	19/26	170	105
059.2008.16	500	24/29	170	105
059.2008.17	500	29/32	170	105
059.2008.18	500	34/35	170	105
059.2008.19	1000	24/29	200	131
059.2008.20	1000	29/32	200	131
059.2008.21	1000	34/35	200	131
059.2008.22	2000	24/29	250	166
059.2008.23	2000	29/32	250	166
059.2008.24	2000	34/35	250	166

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com
<https://www.glacierindia.com>