

Item Name : Flasks Round Bottom

Item Code : 060.2008.01

Item Description :

Two Neck at Angle DIN 12394

Above flasks can also be provided in other joint sizes.

Item Image :



Cat. No.	Nominal capacity ml	Socket size	Side Socket	Approx. overall height mm.
060.2008.01	50	24/29	14/23	115
060.2008.02	50	29/32	14/23	115
060.2008.03	100	24/29	14/23	115
060.2008.04	100	24/29	19/26	115
060.2008.05	100	29/32	14/23	115
060.2008.06	250	24/29	14/23	145
060.2008.07	250	24/29	19/26	145
060.2008.08	250	29/32	14/23	145
060.2008.09	500	24/29	14/23	175
060.2008.10	500	24/29	19/26	175

060.2008.11	500	29/32	14/23	175
060.2008.12	1000	24/29	14/23	210
060.2008.13	1000	24/29	19/26	210
060.2008.14	1000	29/32	14/23	210
060.2008.15	2000	34/35	19/26	223

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>