

Item Name : Flasks Round Bottom Three Neck

Item Code : 061.2008.01

Item Description :



Item Image :

Cat. No.	Nominal capacity ml	Centre Socket	Side Sockets	Approx. overall height mm.
061.2008.01	100	19/26	14/23	115
061.2008.02	100	24/29	14/23	115
061.2008.03	100	24/29	19/26	115
061.2008.03A	100	29/32	14/23	115
061.2008.04	250	19/26	19/26	145
061.2008.05	250	24/29	14/23	145
061.2008.06	250	24/29	19/26	145
061.2008.07	250	29/32	14/23	145
061.2008.08	250	29/32	19/26	145
061.2008.08A	250	29/32	29/32	145
061.2008.09	500	24/29	14/23	175

061.2008.10	500	24/29	19/26	175
061.2008.11	500	29/32	14/23	175
061.2008.12	500	29/32	19/26	175
061.2008.12A	500	29/32	29/32	175
061.2008.13	1000	24/29	24/29	210
061.2008.14	1000	24/29	19/26	210
061.2008.15	1000	29/32	14/23	210
061.2008.16	1000	29/32	19/26	210
061.2008.17	1000	29/32	29/32	210
061.2008.18	2000	24/29	19/26	260
061.2008.19	2000	29/32	24/29	260
061.2008.20	2000	29/32	29/32	290

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com
<https://www.glacierindia.com>