

**Item Name :** Flasks, Round Bottle Four Neck

**Item Code :** 062A.2008.01

**Item Description :**

Centre Neck & three parallel side neck with interchangeable joints. All four necks are vertical. The two side necks are 1800 apart & a third side neck is 900 from other two necks.

**Item Image :**



Cat. No.	Nominal capacity ml	Centre Neck	Side Neck
062A.2008.01	50	24/29	24/29
062A.2008.02	100	24/29	14/23
062A.2008.03	100	24/29	19/26
062A.2008.04	250	24/29	14/23
062A.2008.05	250	24/29	19/26
062A.2008.06	500	24/29	14/23
062A.2008.07	500	24/29	19/26
062A.2008.08	500	24/29	24/29
062A.2008.09	500	34/35	24/29
062A.2008.10	1000	24/29	14/23

062A.2008.11	1000	24/29	19/26
062A.2008.12	1000	24/29	24/29
062A.2008.13	1000	34/35	24/29
062A.2008.14	2000	24/29	14/23
062A.2008.15	2000	24/29	19/26
062A.2008.16	2000	34/35	19/26
062A.2008.17	2000	24/29	24/29
062A.2008.18	5000	34/35	19/26
062A.2008.19	5000	34/35	24/29

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>