

**Item Name :** Flasks Evaporating

**Item Code :** 069.2008.01

**Item Description :**



**Item Image :**

Cat. No.	Nominal capacity ml	Socket	Approx. overall height mm.
069.2008.01	50	29/32	94
069.2008.02	100	24/29	110
110	100	29/32	110
069.2008.04	250	24/29	140
069.2008.05	250	29/32	140
069.2008.06	500	24/29	170
069.2008.07	500	29/32	170
069.2008.08	1000	24/29	210
069.2008.09	1000	29/32	210
069.2008.10	2000	29/32	250

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>