

**Item Name :** 3 Ply Face Mask

**Item Code :** 20007

**Item Description :**

Face Mask 3 ply is made of Polypropylene Spunbond and Meltblown nonwoven fabrics with adjustable nose strip and flat or round elastic ear - loop to give basic protection to the user. It's mainly used in medical, dental, clean room, critical environment, food handling/production and industrial manufacturing.

Quantity : 1 No.

Specifications

Width across face : 175 mm + / - 2 mm

Depth pleated : 95 mm + / - 2 mm

Depth in full : 175 mm + / - 2 mm

Nose strip : 95 mm + / - 2 mm

Loop size : 155 mm + / - 5 mm

Components

1st layer : Polypropylene Spunbond nonwoven

2nd layer : Meltblown Nonwoven

3rd layer : Polypropylene Spunbond nonwoven

Nose strip : PE wire

Ear-Loop : Flat elastic or round elastic (latex free)

Benefits and Standards

1. Excellent breath ability.
2. Non irritating, fiber glass free fabric.
3. Hypoallergenic .
4. Non-woven spunbond 100% Polypropylene.
5. Splash resistant
6. Hygienic

**Item Image :**



Cat. No.

20007

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>