

**Item Name :** Bouffant Cap

**Item Code :** 20009

**Item Description :**

Features of Disposable Non-Woven Cap :

Adjustable soft thread elastic for comfort fittings, avoids escape of hair, dandruff and micro organisms during process.

Technical Specification of Disposable Non-Woven Cap :

These caps are made from spun bond polymer non-woven fabric. This material is preferred for its light weight. Its gravity is lighter than hydrophobic PP fibers so it has excellent water repellency and air permeability resulting in providing comfort & dryness resulting in anti bacterial and anti fungal resistance.

Quantity Material Colour Size Thickness

1 no. PP Spun Blue, Green Free Size 30-40 GSM

Bond Fabric White

**Item Image :**



Cat. No.

20009