


Item Name : Test Tubes Round Bottom with rim, Plain Boro 3.3

Item Code : 097.2008.01

Item Description :



Item Image :

Cat. No.	Capacity ml	O.D. (mm) x Height	Wall Thickness
097.2008.01	3	10x75	1.0
097.2008.02	5	12x75	1.0
097.2008.03	8	12x100	1.0
097.2008.04	10	15x125	1.2
097.2008.05	25	18x150	1.2
097.2008.06	50	25x150	1.2
097.2008.07	75	25x200	1.2
097.2008.08	100	32 x 200	1.8

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>