

**Item Name :** Test Tubes Round Bottom without Rim

**Item Code :** 098.2008.01

**Item Description :**



**Item Image :**

Cat. No.	Capacity ml	O.D. (mm) x Height	Wall Thickness
098.2008.01	3	10x75	1.0
098.2008.02	5	12x75	1.0
098.2008.03	8	12x100	1.0
098.2008.04	10	15x125	1.2
098.2008.05	25	18x150	1.2
098.2008.06	50	25x150	1.2
098.2008.07	75	25x200	1.2
098.2008.08	100	32 x 200	1.8

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

[info@glacierindia.com](mailto:info@glacierindia.com)

<https://www.glacierindia.com>