

Item Name : Student Microscope With Fixed Condenser

Item Code : GW - 1105

Item Description :

Model	GW - 1105
All Make Optics Having Features	Para focal, precentred, high numerical aperture, high resolution
Special Coating Having No Effect Of	Moisture, fungus, fogginess or any other weather effect with passage of time.
Magnification	50x to 675x
Stand	All metal constructed body with properly balanced , inclinable arm up to 90 degree for easy viewing
Focussing	Separate coarse fine focusing knobs
Stage	Full square stage 112 x 112 mm size with stage clips
Illumination	Plano concave mirror on gimble mount
Nose piece	Triple revolving
Objectives	10x & 40/45x (s. L.) (all achromatic).
Eye pieces	10x & 15x
Condenser	Fixed abbe type condenser with iris diaphragm
Packing	Nicely packed in box.

Item Image :



Model

GW - 1105

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>