

1141.2008.01 Round Screw

1141.2008.02 Thumb Screw

GLACIER INDIA

 **Glacier**
INDIA
www.glacierindia.com

(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)



Item Name : Boss Head (Cast Iron)

Item Code : 1141.2008.03

Item Description :

For use with retort and condenser clamp. painted, with strong screw to hold up to 20 mm light. Different types of screw available. As per DIN 12895.

Item Image :



Cat. No. Type

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>