

**Item Name :** Student Microscope with Movable Condenser

**Item Code :** GW - 1106

**Item Description :**

<b>Model</b>	GW - 1106
<b>All Optics Having Features</b>	Parafoval, precentred, high numerical aperture, high resolution
<b>Special Coating Having No Effect Of</b>	Moisture, fungus, fogginess or any other weather effect with passage of time.
<b>Magnification</b>	50x to 675x
<b>Stand</b>	All metal constructed body with properly balanced, inclinable arm upto 90 degree for easy viewing
<b>Focussing</b>	Separate coarse fine focusing knobs
<b>Stage</b>	Full square stage 112x 112mm size with stage clips
<b>Illumination</b>	Planoconcave mirror on gimble mount
<b>Nose piece</b>	Triplerevolving
<b>Objectives</b>	10x & 40/45x (s. L. ) (all achromatic).
<b>Eye pieces</b>	10x & 15x
<b>Condenser</b>	Fixed abbe type condenser with iris diaphragm
<b>Packing</b>	Nicely packed in box.

**Item Image :**



Model

GW - 1106

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

[info@glacierindia.com](mailto:info@glacierindia.com)

<https://www.glacierindia.com>