

Item Name : Student Medical Microscope For Educational Purpose

Item Code : GW-1108

Item Description :

Model	GW-1108
Magnification	20x to 1500x
Stand	Robust and stable, inclinable through 90 degrees. Mechanical tube length 160mm.
Focussing	Separate knobs for coarse and fine motion.
Stage	Fixed square stage 120 x 120 mm size with mechanical stage having travel 30 x 76mm with graduated scales and vernier.
Illumination	Plano concave mirror on gimble mount. Optionally available attachable 20W for 220V or 6V, 20W with transformer.
Nose piece	Quadruple-revolving
Objectives	10x & 45x (s. L.) with 100x oil immersion.
Eye pieces	Huygenian 5x or 15x, wf 10x.
Condenser	Abben. A. 1. 25 with iris diaphragm focusing up and down by rack & pinion.
Packing	Nicely packed in box.

Item Image :



Model

GW-1108



44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>