


**Item Name :** Student Monocular Microscope for Educational Purpose

**Item Code :** GW-1110

**Item Description :**

|                      |  |
|----------------------|--|
| <b>Model</b>         | GW-1110  |
| <b>Observation</b>   | Monocular- interchangeable with binocular or vertical photo tube.  |
| <b>Magnification</b> | 20x to 1500  |
| <b>Stand</b>         | One piece die cast arm with ball bearing guide ways providing extra smooth and frictionless motion to the fine focusing mechanism without any back lash. |
| <b>Focussing</b>     | Coarse and fine focusing by conveniently placed separate knobs. 1 div. = 0.002 mm.   |
| <b>Stage</b>         | Mechanical stage of standard stage size 120x 140 mm.   |
| <b>Illumination</b>  | Plano concave mirror available with optional lamp for 6v/20w or direct 220v attachable in place of the mirror mount. .                                   |
| <b>Nose piece</b>    | Quadruple- revolving on smooth bearings.   |
| <b>Objectives</b>    | Achromat 4x,10x, 40x (s. L. ) & 100x oil immersion.  |
| <b>Eye pieces</b>    | Huygenian 5x or 15x, w/ 10x.   |
| <b>Condenser</b>     | Bright field, abbe's n. A. 1. 25 focusing up and down by rack & pinion.  |
| <b>Packing</b>       | Nicely packed in box.  |

**Item Image :**



Model

GW-1110



# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>