


Item Name : Student Monocular Microscope for Educational Purpose

Item Code : GW-1110

Item Description :

Model	GW-1110
Observation	Monocular- interchangeable with binocular or vertical photo tube.
Magnification	20x to 1500
Stand	One piece die cast arm with ball bearing guide ways providing extra smooth and frictionless motion to the fine focusing mechanism without any back lash.
Focussing	Coarse and fine focusing by conveniently placed separate knobs. 1 div. = 0.002 mm.
Stage	Mechanical stage of standard stage size 120x 140 mm.
Illumination	Plano concave mirror available with optional lamp for 6v/20w or direct 220v attachable in place of the mirror mount. .
Nose piece	Quadruple- revolving on smooth bearings.
Objectives	Achromat 4x,10x, 40x (s. L.) & 100x oil immersion.
Eye pieces	Huygenian 5x or 15x, w/ 10x.
Condenser	Bright field, abbe's n. A. 1. 25 focusing up and down by rack & pinion.
Packing	Nicely packed in box.

Image not found or type unknown



Item Image :

Model

GW-1110

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>