

**Item Name :** Pathological Monocular Microscope

**Item Code :** GW-1111

**Item Description :**

<b>Model no.</b>	GW-1111
<b>Observationhead</b>	The monocular version with 45 degree inclined monocular tube rotate-able through 360 degree.
<b>Mainbody</b>	All part pressure die cast to provide long life and good mechanical strength.
<b>Eye pieces</b>	Wide field 10x or 15x.
<b>Objectives</b>	4x, 10x, 40/45x (s. L.), 100x oil immersion (s. L.)
<b>Nose pieces</b>	Quadruple nose piece revolves on ball bearing.
<b>Mechanical stage</b>	Co-axial low positioned mechanical x-y stage of size 120x 125mm with stage fine vernier graduation.
<b>Condenser</b>	Abbe condenser n. A. 1. 25 having a spherical lens, adjustable iris diaphragm.
<b>Illumination</b>	Built in base illumination system with 6v-20w halogen bulb and plano concave mirror attachment suitable for 220v, 50hz ac supply
<b>Packing</b>	Nicely packed in box.

**Item Image :**



Model no.

GW-1111

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

[info@glacierindia.com](mailto:info@glacierindia.com)

<https://www.glacierindia.com>