

Item Name : Binocular Pathological Microscope

Item Code : GW-1112

Item Description :

Model no.	GW-1112
Observationhead	Binocular version with 45 degree inclined monocular tube rotate-able through 360 degree. Interpupillary distance adjustment ring 54-74 mm& dioptic adjustment
Mainbody	All part pressure diecast to provide long life and good mechanical strength.
Eyeieces	Widefield 10x or 15x.
Objectives	4x,10x,40/45x (S. L.), 100x oil immersion(S. L.)
Nosepieces	Quadruple nosepiece revolves on ball bearing.
Mechanical stage	Co-axial low positioned mechanical x-y stage of size 120x 125mm with stage fine vernier graduation.
Condenser	Abbe condenser n. A. 1. 25 having a spherical lens, adjustable iris diaphragm.
Illumination	Built in base illumination system with 6v-20w halogen bulb and plano concave mirror attachment suitable for 220v, 50hz ac supply
Packing	Nicely packed in box.

Item Image :



Model no.

GW-1112

A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>