

**Item Name :** Multi Viewing Microscope

**Item Code :** GW-1150

**Item Description :**

<b>Model</b>	GW-1150
<b>Optical system</b>	Infinite optical system
<b>Viewing head</b>	Compensation free trinocular head, inclined at 30 °, 360 ° rotatable; interpupillary distance Compensation free binocular head, inclined at 30 °, 360 ° rotatable; interpupillary distance 48-75 mm.
<b>Eye piece</b>	Extra wide field eyepiece wf 10x/20.
<b>Nosepiece</b>	Inverse quintuple nosepiece
<b>Objective</b>	Infinite plan achromatic objectives 4x, 10x, 40x, 100x.
<b>Condenser</b>	Swing condenser n. A. 0. 9/0. 25
<b>Focusing</b>	Coaxial coarse & fine focus adjustment, fine division 0. 002 mm.
<b>Stage</b>	Double layers mechanical stage 185 x 142 mm, moving range 75 x 55 mm.
<b>Kohler illumination</b>	External illumination, halogen lamp 12v 50w.
<b>Pointer</b>	Green led pointer, brightness adjustable
<b>Video adapter</b>	C mount 1x C mount 0. 5x

Item Image :



Model

GW-1150

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>