

Item Name : Multi Viewing Microscope

Item Code : GW-1150

Item Description :

Model	GW-1150
Optical system	Infinite optical system
Viewing head	Compensation free trinocular head, inclined at 30 °, 360 ° rotatable; interpupillary distance Compensation free binocular head, inclined at 30 °, 360 ° rotatable; interpupillary distance 48-75 mm.
Eye piece	Extra wide field eyepiece wf 10x/20.
Nosepiece	Inverse quintuple nosepiece
Objective	Infinite plan achromatic objectives 4x, 10x, 40x, 100x.
Condenser	Swing condenser n. A. 0. 9/0. 25
Focusing	Coaxial coarse & fine focus adjustment, fine division 0. 002 mm.
Stage	Double layers mechanical stage 185 x 142 mm, moving range 75 x 55 mm.
Kohler illumination	External illumination, halogen lamp 12v 50w.
Pointer	Green led pointer, brightness adjustable
Video adapter	C mount 1x C mount 0. 5x

Item Image :



Model

GW-1150

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>