

Cat No.

2063.2008.01



GLACIER INDIA
(ISO,CE,WHO,GMP,FDA 17025:2017 Certified)



Item Name : Filter Paper With Different Sizes

Item Code : 2063.2008.01

Item Description :

- 5.5cm Diameter, Pore Size: 20-25 μ m
- Grade 4 Qualitative Filter Paper
- Used in qualitative analytical techniques to determine and identify materials
- Extremely fast filtering; provides excellent retention of coarse particles (20-25 μ m) and gelatinous precipitates such as ferric hydroxide and aluminum hydroxide
- Used when high flow rates in air pollution monitoring are required and the collection of fine particles is not critical • Available in different sizes

Item Image :



44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>