

Item Name : Research Trinocular Ore Microscope

Item Code : GW - Research ore

Item Description :

Model	GW - Research ore
Observation head	Flat trinocular head inclined at 45°
Eyepiece	Eyepiece 10x20mm Cross partition eyepiece 10x20mm Grid eyepiece 10x20mm
Objectives	4x, 10x, 25/40x, 60/63x
Nosepiece	Quadruple nosepiece
Mechanical stage	Revolving round stage diameter 160mm (approx. 360° rotatable & graduated in 1°) Least count 0. 1° with vernier Compensator-gypsum 1°, mica ¼°, quartz wedge
Illumination	12v/20w lamp brightness adjustable Power supply 110v-220v.
Optional accessories	1) Digital high resolution photographic camera with microscope adopter.

Item Image :



Model

GW - Research ore

44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>