

**Item Name :** Toolmaker Microscope

**Item Code :** GW-1101

**Item Description :**

For making various types of measurements of angles , Thickness radius, diameter etc.

Magnification; 30x with objectives 2x & Eyepiece 15x with field of view 12mm & working distance about 100mm.

Microscope tube inclined at 45 degree. Measuring stage 150x150mm with travels in X-Y direction. Travel range 50 x 50mm having micrometer

Rotatory stage graduation 0-360 degree with vernier 6 minutes complete provided with inclined Transmitted light illumination having separate

Optional Accessories Rotatory measurement stage: Circular stage rotatable through 0-360 degree having vernier

**Item Image :**



Model

GW-1101