

**Item Name :** Constant Deviation Spectrograph

**Item Code :** GW-11050

**Item Description :**

**Base:** Heavy stable base having a pillar on which mounted a platform collimator/telescope prism housing. **Prism**

**Telescope:** Fitted on the platform, focussed by rack/pinion, achromatic objective with light resolution, eyepiece

**Collimator:** With achromatic objective having high resolution operateable by rack & pinion, having slit variable

**Prism:**

Glass prism e. D. F. Glass prism d. E. D. F.

Extension mount for keeping I. G. Plate/edser butle Plate/ fabry perot etalon.

**Extra attachments for telescope.**

- Attachment for converting c. D. S. Into monochromator.
- Micrometer eyepiece – filer type reading 0. 001 cms
- Camera for above with plate holder.
- Plate holder for camera.

**Item Image :**



Model

GW-11050

# 44, Vikas puri, P.O. Indraprastha, Jagadhri Road, @glacierindia.com  
A/Cantt-133 006. Haryana, India. <https://www.glacierindia.com>