

**Item Name :** Research Inverted Metallurgical Microscope

**Item Code :** GW - 2001

**Item Description :**

<b>GW - 2001</b>	
<b>Observation head</b>	Side & top with eyepiece tube Binocular observation head inclined at 30°, interpupillary distances adjustment range from 55 to 75mm.
<b>Main body</b>	All parts die cast to provide long life & good mechanical strength. Co-axial focusing control knobs, tension adjustment ring & pre focusing lock.
<b>Eyepiece</b>	10x (16mm).
<b>Objectives</b>	PL10x, PL25x, PL40x, PL60x (optional).
<b>Nosepiece</b>	Quadruple nosepiece revolves on ball bearing.
<b>Mechanical stage</b>	Co-axial mechanical stage size 150x200mm cross travel range 15x15mm.
<b>Illumination</b>	Koehler illumination system, 6V20W halogen lamp with adjustable brightness.
<b>Packing</b>	Nicely packed in box.

**Item Image :**



Model

GW - 2001

# 44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,  
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>