

Item Name : Research Inverted Metallurgical Microscope

Item Code : GW - 2001

Item Description :

GW - 2001	
Observation head	Side & top with eyepiece tube Binocular observation head inclined at 30°, interpupillary distances adjustment range from 55 to 75mm.
Main body	All parts die cast to provide long life & good mechanical strength. Co-axial focusing control knobs, tension adjustment ring & pre focusing lock.
Eyepiece	10x (16mm).
Objectives	PL10x, PL25x, PL40x, PL60x (optional).
Nosepiece	Quadruple nosepiece revolves on ball bearing.
Mechanical stage	Co-axial mechanical stage size 150x200mm cross travel range 15x15mm.
Illumination	Koehler illumination system, 6V20W halogen lamp with adjustable brightness.
Packing	Nicely packed in box.

Item Image :



Model

GW - 2001



44, Vikas puri, P.O.Ind.Estate, Jagadhri Road,
A/Cantt-133 006. Haryana, India.

+0171 – 2698107

info@glacierindia.com

<https://www.glacierindia.com>