

**Item Name :** Ore Polarising Microscope

**Item Code :** GW-120OR

**Item Description :**

Microscope frame for transmission and reflection mode polarization microscopy with built-in filter holders

Infinity Corrected

Plan achromatic, polarizing objectives with magnifications 4X POL N.A. 0.13 W.D. 17.0 mm

10X POL N.A. 0.3 W.D. 10.0 mm

20X POL N.A. 0.5 W.D. 2.1 mm

40X POL N.A. 0.75 W.D. 1.1 mm

Wide field 10X eyepiece (02 numbers), with double cross line focusable . Minimum 22 mm

Fixed Trinocular Tube.

High Power LED Light with power supply .

High Power LED light with power supply Built in Transmitted light filters. Minimum 4 position centrable & detachable nosepiece .

Stroke: Min 22mm. Fine Stroke: Min 100µm Minimum

**Item Image :**



Model

GW-120OR